

Avalanche®-Omni Transfection Reagent

Cat. No. EZT-OMNI-1 Size: 0.75 ml Store at 4°C

1.5 ml

Description

Avalanche®-Omni Transfection Reagent is a proprietary lipid-polymer mixture that forms lipopolyplexes with nucleic acids. It is an exceptionally powerful and versatile next-generation DNA and siRNA transfection reagent for day-to-day experiments. It has been shown to effectively transfect the broadest spectrum of adherent and suspension cells. No other transfection reagents can match the efficiency, broad spectrum capabilities, convenience, and gentleness of Avalanche®-Omni Transfection Reagent for transfection of primary cells and other hard-to-transfect cells. Because of the high transfection efficiency, much lesser amounts of DNA and Avalanche®-Omni are needed to achieve high transfection efficiency as compared to most other transfection reagents. It is best to use the smallest amounts of DNA and transfection reagent possible in order to maintain good cell growth and viability, and generate physiologically relevant data you can trust.

Features:

- Exceptional transfection efficiency in the broadest range of cell types including difficult-to-transfect cell lines and primary cells with super high levels of gene expression or knockdown.
- Economical: High efficiency means smaller amounts of the nucleic acid and reagent are needed. 1.5 ml of Avalanche®-Omni Reagent is sufficient for about 3,500 transfections in 24-well plates, or about 700 transfections using 6-well plates!
- Extremely gentle to cells: The chemical features of our proprietary formula and the smaller amounts needed for transfection ensure very low toxicity levels.
- Simplest procedure: simply mix with nucleic acid and add to cell culture.
- Superior performance for co-transfection of siRNA and plasmid DNA
- Proven efficacy in the presence of serum eliminates the need to change media following transfection
- 100% animal origin- free
- Reliable performance for high-throughput applications
- The best choice for establishing stable cell lines

BEFORE YOU START:

Important Tips for Optimal Transfection

- 1. Prepare high-quality plasmid DNA at $0.5-5~\mu g/\mu l$ in deionized water or TE buffer. Make sure the plasmids are endotoxin-free and have A260/280 absorbance ratio of 1.8-2.0. A GFP (green fluorescent protein) plasmid can be used to determine transfection efficiency.
- 2. Use Opti-MEM® I Reduced Serum Medium (Life Technologies) or regular DMEM without serum to make Avalanche®-Omni and nucleic acid mix (Only Opti-MEM® I will be mentioned in the remaining parts of the protocol for simplification purposes). Do not use NaCl₂ solution or PBS.
- 3. Maintain the same seeding conditions between experiments. Use low-passage cells; make sure that cells are healthy and greater than 90% viable before transfection.
- 4. It is important to have the cells in proliferation state and 70-90% confluence at the time of DNA transfection.
- 5. Avalanche®-Omni Transfection Reagent is extremely gentle to cells. However, transfection process will impose stress on cells, no matter what type of transfection methods you use. The trick is to get the balance between transfection efficiency and cell viability. Increasing the number of cells plated per well or decreasing Avalanche®-Omni and DNA amount will minimize the effect of transfection on cell growth and viability. With careful optimization, as described in page 3 and 4, this can be achieved while keeping the highest transfection efficiency.
- 6. Don't use antibiotics in the culture medium during the first 24 hours of transfection.

Protocols

1 DNA Transfection

1.1 Cell Seeding

For optimal DNA transfection conditions, we recommend using cells which are 70% to 90% confluent at the time of transfection. Typically, for experiments in 24-well plates, 50,000-80,000 adherent cells are seeded per well in 0.5 ml of cell growth medium **without antibiotics** 24 h prior to transfection. For the different culture formats, refer to Table 1.

2

Table 1. Recommended number of cells to seed the day before transfection in culture medium without antibiotics

Culture vessel	Number of Adherent cells to seed (Suspension Cells)	Surface area per well (cm²)	Volume of medium per well to seed the cells (ml)	
96-well	7,500-10,000 (4x10 ⁴)	0.3	0.1	
24-well	50,000-80,000 (2x10 ⁵)	1.9	0.5	
12-well	80,000-150,000 (4x10 ⁵)	3.8	1	
6-well/35 mm	150,000-250,000 (8x10 ⁵⁾	9.4	2	
60 mm/flask 25 cm ²	250,000-800,000 (2x10 ⁶)	25-28	5	
100 mm/flask 75 cm ²	1x10 ⁶ -2x10 ⁶ (6x10 ⁶)	75-78.5	10	
150 mm/flask 175 cm ² 2x10 ⁶ -5x10 ⁶ (1.3x10 ⁷)		153-175	25	

1.2 DNA Transfection


If this is the first time that you are using Avalanche®-Omni on a specific type of cells, follow the following procedures and Figure 1 for optimization (The optimization procedures are extremely important for successful transfection. Since different types of cells have different sensitivity to Avalanche®-Omni, the amount of Avalanche®-Omni as well as the amount of DNA needed for maximum transfection on different types of cells may differ dramatically).

As an example, the following procedures and Figure 1 are for optimization on 24-well plate.

3

- 1. Bring Avalanche®-Omni and serumfree medium (OptiMEM®I) to room temperature before starting.
- 2. Add 200μl of OptiMEM®I into two 1.5ml tubes (Tube A and B).
- 3. Add 4.0 μg of DNA to Tube A and 8.0 μg of DNA to Tube B.
- 4. Mix Avalanche®-Omni prior to preparing complexes. Add 6.4 μl of Avalanche® -Omni to both Tube A and Tube B containing 200μl of different concentration of DNA solution (20.0 μg/ml and 40.0 μg/ml respectively). Vortex for 2-3 seconds. Incubate each tube for 15 minutes at room temperature.
- 5. After incubation, add 3.12, 6.25, 12.5, 25, 50, and 100μl of DNA/Avalanche®-Omni complexes dropwise directly to the corresponding wells of the 24-well cell culture plate (See Figure 1). Swirl plate gently.
- 6. Incubate the cells at 37°C in a CO₂ incubator
- 7. Change growth medium 12-16 hours later.
- 8. Expression of reporter gene activity should generally be assessed at 20-48 hours post-transfection. GFP expression is maximal at 40-48 hours post-transfection.

Figure 1.



After you have completed the optimization steps, choose the amount of DNA and Avalanche®-Omni that gave you the optimal balance of potency & low cytotoxicity (which usually is the lowest dose that gave you the same high transfection efficiency as other higher doses did) for all of your future experiments on this specific cell type.

1.3 Scale Up or Down Transfections

Use Table 2 to scale the amount of DNA/Avalanche®-Omni for your transfection experiment.

Table 2. Scaling Up or Down Transfection Instruction

Culture Vessel	Multiplication factor ¹		
96-well	0.17		
48-well	0.50		
24-well	1.00		
12-well	2.00		
6-well	5.00		
60-mm	11.05		
10-cm	28.95		
T75	39.47		

¹After determining the optimum amount of DNA/Avalanche®-Omni from the optimization procedures on the above 24-well plate, use the multiplication factor to determine the DNA and Avalanche®-Omni amount needed for your new plate format.

2 siRNA Transfection

2.1 Cell Seeding

For optimal siRNA transfection conditions, we recommend using cells which are 50% confluent at the time of transfection. Typically, for experiments in 6-well plates, 100,000 to 150,000 cells are seeded per well in 2 ml of growth medium without antibiotics 24 h prior to transfection. For other culture formats, refer to Table 3.

Table 3. Recommended number of cells to seed the day before transfection in culture medium without antibiotics

Culture vessel	Number of adherent cells to seed	Surface area per well (cm²)	Medium per well to seed the cells (ml)	
24-well	25,000-40,000	1.9	0.5	
12-well	50,000-80,000	3.8	1	
6-well/35 mm	100,000-150,000	9.4	2	
60 mm/flask 25 cm ²	200,000-500,000	25-28	5	
100 mm/flask 75 cm ²	0.5x10 ⁶ -1x10 ⁶	75-78.5	10	

5

2.2 siRNA Transfection

For optimal siRNA-mediated silencing, we recommend using 10 to 50 nM siRNA (final concentration). The following conditions are given per well of a 6-well plate. For other culture formats, please refer to Table 4.

- 1. Dilute 22 to 110 pmoles siRNA (final concentration: 10 to 50 nM) into 200 μ l of Opti-MEM®I. Mix by vortexing.
- 2. Briefly vortex Avalanche®-Omni, and add 1.0-5.0 μ l into the diluted siRNA. Immediately vortex for 10 s.
- 3. Incubate for 15 min at RT.
- 4. Add the transfection mixture drop-wise into each well.
- 5. Gently rock the plates back and forth and from side to side, and return the plate to the 37°C CO2 incubator.
- 6. Analyze after incubating for 24 h or longer.

Table 4. siRNA transfection guidelines according to the cell culture vessel per well

Culture Vessel	siRNA (pmole) 10 nM	siRNA (pmole) 50 nM	Avalanche®- Omni (μl)	Opti-MEM or DMEM (µl))	Growth medium (ml)	Final Volume in the well (ml)
24-well	5.5	27.5	*0.2-1.0	50	0.5	0.55
12-well	11	55	*0.4-2.0	100	1	1.1
6-well/ 35 mm	22	110	1.0-5.0	200	2	2.2
60 mm/ flask 25 cm ²	44	220	2.3-11.5	400	4	4.4
100 mm/ flask 75 cm ²	121	605	5.8-29	1100	11	12.1

^{*} Dilute Avalanche®-Omni 1:5 with H_2O prior to application (4 μ l reagent + 16 μ l H_2O), and then use 5 times of the volume in the table for accurate pipetting.



Avalanche®-Omni Transfection Reagent

Intended Use:

All Avalanche® Series Transfection Reagents are for research use only, not intended for any animal or human therapeutic or diagnostic use.

Avalanche® is the trademark of EZ Biosystems LLC. All rights reserved.

7

Avalanche® series Transfection Reagents are sold to the Buyer with a limited license for Research Use Only and is not for clinical, therapeutic or diagnostic use in humans or animals.

Avalanche® series Transfection Reagents may not be re-packaged or re-sold without written permission from EZ Biosystems LLC. The buyer agrees not to attempt to reverse engineer, reconstruct, synthesize or otherwise modify EZ Biosystems Avalanche® products. A license from EZ Biosystems LLC is required for commercial application of Avalanche® series Transfection Reagents. For obtaining a license to use those products for commercial application, contact EZ Biosystems LLC via business@ezbiosystems.com. For full terms and conditions, visit: www.ezbiosystems.com.